

40%

Of women in the U.S.
have dense breasts

40%

Of cancers are missed
by mammography in
extremely dense
breasts

30%

Of cancer missed
by mammography
in **heterogeneously**
dense breasts

4-6 Xs

Greater risk of developing breast
cancer for women with dense
breasts compared to
women with fatty breasts.

BIRADS c & d

Cancer
Detection ↓ ↑ Cancer
Risk

Your patients deserve an **efficient** and **effective** solution for
dense breast screening.

Talk to them about **SoftVue™ 3D Whole Breast
Ultrasound Tomography.**



844-SOFTVUE • DelphinusMT.com



The Technology

- **Patented circular-array transducer** and **Sequir™ Breast Interface** provide the basis for the image acquisition.
- Three sound signals: **Reflection, Sound Speed** and **Attenuation** captured during image acquisition providing tissue structure and characterization properties.
- Reconstructed into **four image sequences** for interpretation in combination with mammography and presented in the coronal plane allowing radiologists to scroll through the entire volume of the breast.



The Scan



- **No Compression** or radiation
- **Comfortable** bed & warm water coupling agent
- Performed on **same day as mammography**
- **2-3 min average** scan time per breast
- **95%** of women would recommend

The Clinical Evidence

The FDA approved **SoftVue™** 3D Whole Breast Ultrasound Tomography System for use as an **adjunct to digital mammography in the screening of asymptomatic women with dense breast tissue**. Clinical evidence has demonstrated that SoftVue enhances dense breast screening identifying more cancers with greater accuracy and potentially fewer biopsies than full field digital mammography alone.

20% Increase in Sensitivity

8% Increase in Specificity

1. Berg, W.A. Breast MRI for "the Masses". Eur Radiol 32, 4034–4035 (2022);

2. Densebreast-info, inc.. DenseBreast. (2022, May 18). Retrieved October 7, 2022, from <https://densebreast-info.org/providers-faqs/do-dense-breasts-affect-the-risk-of-developing-breast-cancer/>

3. Delphinus Medical Technologies, Inc (2021). SOFTVUE™ 3D WHOLE BREAST ULTRASOUND TOMOGRAPHY SYSTEM (SOFTVUE™) PHYSICIAN LABELING. https://www.accessdata.fda.gov/cdrh_docs/pdf20/P200040C.pdf